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👉 **Lưu ý:**

- Before performing the isolation protocol, **carefully** read the **note** in protocol
- **Follow** the **compliance** when **pipetting**

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### Wizard® Genomic DNA Purification Kit

INSTRUCTIONS FOR USE OF PRODUCTS A1120, A1123, A1125 AND A1620.

Quick  
PROTOCOL

#### Isolation of Genomic DNA from Whole Blood

Sample Size	Lysis Solution		Protein Precipitation Solution	Isopropanol	DNA Rehydration Solution
	Cell	Nuclei			
300µl	900µl	300µl	100µl	300µl	100µl
1ml	3ml	1ml	330µl	1ml	150µl
3ml	9ml	3ml	1ml	3ml	250µl
10ml	30ml	10ml	3.3ml	10ml	800µl

As little as 20µl can be processed with this system. Please see Technical Manual #TM050, Section 3.C.

#### Red Blood Cell Lysis

- Using volumes from the table above, combine the appropriate volumes of Cell Lysis Solution and blood. Mix by inversion.
- Incubate for 10 minutes at room temperature.
- Centrifuge:  
≤300µl sample 13,000–16,000 × g\*; 20 seconds  
1–10ml sample 2,000 × g; 10 minutes
- Discard supernatant. Vortex pellet.


#### Nuclei Lysis and Protein Precipitation

- Using volumes from the table above, add Nuclei Lysis Solution and mix by inversion.
- Add Protein Precipitation Solution; vortex for 20 seconds.
- Centrifuge:  
≤300µl sample 13,000–16,000 × g\*; 3 minutes  
1–10ml sample 2,000 × g; 10 minutes

#### DNA Precipitation and Rehydration

- Transfer supernatant to a new tube containing isopropanol (using volumes from table above). Mix.
- Centrifuge:  
≤300µl sample 13,000–16,000 × g\*; 1 minute  
1–10ml sample 2,000 × g; 1 minute
- Discard supernatant. Add 70% ethanol (same volume as isopropanol).
- Centrifuge as in Step 9.
- Aspirate the ethanol and air-dry the pellet (10–15 minutes).
- Rehydrate the DNA in the appropriate volume of DNA Rehydration Solution for 1 hour at 65°C or overnight at 4°C.

\*Maximum speed on a microcentrifuge.

2818MAT1\_JA  
2/18/04

### Wizard® Genomic DNA Purification Kit

INSTRUCTIONS FOR USE OF PRODUCTS A1120, A1123, A1125 AND A1620.

Quick  
PROTOCOL

#### Isolation of Genomic DNA from Animal Tissue and Tissue Culture Cells

##### Prepare Tissues

**Tissue Culture Cells:** Centrifuge at 13,000–16,000 × g\* for 10 seconds. Wash the cell pellet with PBS, vortex and then add 600µl of Nuclei Lysis Solution and mix by pipetting.

**Animal Tissue:** Add 10–20mg of fresh or thawed tissue to 600µl of chilled Nuclei Lysis Solution and homogenize for 10 seconds. Alternatively, use 10–20mg of ground tissue. Incubate at 65°C for 15–30 minutes.

**Mouse Tail:** Add 600µl of chilled EDTA/Nuclei Lysis Solution to 0.5–1cm of fresh or thawed mouse tail. Add 17.5µl of 20mg/ml Proteinase K and incubate overnight at 55°C with gentle shaking.

##### Lysis and Protein Precipitation


- Add 3µl of RNase Solution to the cell or animal tissue nuclei lysate and mix. Incubate for 15–30 minutes at 37°C. Cool to room temperature.
- Add 200µl of Protein Precipitation Solution. Vortex and chill on ice for 5 minutes.
- Centrifuge at 13,000–16,000 × g\* for 4 minutes.

##### DNA Precipitation and Rehydration

- Transfer supernatant to a fresh tube containing 600µl of room temperature isopropanol.
- Mix gently by inversion.
- Centrifuge at 13,000–16,000 × g\* for 1 minute.
- Remove supernatant and add 600µl of room temperature 70% ethanol. Mix.
- Centrifuge as in Step 6.
- Aspirate the ethanol and air-dry the pellet for 15 minutes.
- Rehydrate the DNA in 100µl of DNA Rehydration Solution for 1 hour at 65°C or overnight at 4°C.

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