

BAM Appendix 2: Most Probable Number from Serial Dilutions

October 2010

Bacteriological Analytical Manual

Appendix 2 Most Probable Number from Serial Dilutions

Author: Robert Blodgett (retired)

For additional Information, contact: [Stuart Chirtel \(mailto:stuart.chirtel\)](mailto:stuart.chirtel) or [Guodong Zhang \(mailto:guodong.zhang@fda.hhs.gov\)](mailto:guodong.zhang@fda.hhs.gov)

Revision History:

- April 2015: Contact for this Appendix was updated
- October 2010: Equation for most probable number (MPN) replaced with graphical version; Added explanatory note for the downloadable spreadsheet.
- July 2003: Added Table 5. for 10 tubes at 10 ml inocula and link to spreadsheet.

Background

References

Tables

Background

Serial dilution tests measure the concentration of a target microbe in a sample with an estimate called the most probable number (MPN). The MPN is particularly useful for low concentrations of organisms (<100/g), especially in milk and water, and for those foods whose particulate matter may interfere with accurate colony counts. The following background observations are adapted and extended from the article on MPN by James T. Peeler and Foster D. McClure in the *Bacteriological Analytical Manual* (BAM), 7th edition.

Only viable organisms are enumerated by the MPN determination. If, in the microbiologist's experience, the bacteria in the prepared sample in question can be found attached in chains that are not separated by the preparation and dilution, the MPN should be judged as an estimate of growth units (GUs) or colony-forming units (CFUs) instead of individual bacteria. For simplicity, however, this appendix will speak of these GUs or CFUs as individual bacteria. If a confirmation test involves selecting colonies to test, then a statistical adjustment not discussed in this appendix should be used (see Blodgett 2005a.)

The following assumptions are necessary to support the MPN method. The bacteria are distributed randomly within the sample. The bacteria are separate, not clustered together, and they do not repel each other. Every tube (or plate, etc.) whose inoculum contains even one viable organism will produce detectable growth or change. The individual tubes of the sample are independent.

The essence of the MPN method is to dilute the sample to such a degree that inocula in the tubes will sometimes but not always contain viable organisms. The "outcome", i.e., the number of tubes and the number of tubes with growth at each dilution, will imply an estimate of the original, undiluted concentration of bacteria in the sample. In order to obtain estimates over a broad range of possible concentrations, microbiologists use serial dilutions incubating tubes at several dilutions.

The MPN is the number which makes the observed outcome most probable. It is the solution for λ , concentration, in the following equation

$$\sum_{j=1}^k \frac{g_j m_j}{1 - \exp(-\lambda m_j)} = \sum_{j=1}^k t_j m_j$$

where $\exp(x)$ means e^x , and

K denotes the number of dilutions,

g_j denotes the number of positive (or growth) tubes in the j th dilution,

m_j denotes the amount of the original sample put in each tube in the j th dilution,

t_j denotes the number of tubes in the j th dilution.

In general, this equation can be solved by iteration.

McCrary (1915) published the first accurate estimation of the number of viable bacteria by the MPN method. Halvorson and Ziegler (1933), Eisenhart and Wilson (1943), and Cochran (1950) published articles on the statistical foundations of the MPN. Woodward (1957) recommended that MPN tables should omit those combinations of positive tubes (high for low concentrations and low for high concentrations) that are so improbable that they raise concerns about laboratory error or contamination. De Man (1983) published a confidence interval method that was modified to make the tables for this appendix.

Confidence Intervals

The 95 percent confidence intervals in the tables have the following meaning:

Before the tubes are inoculated, the chance is at least 95 percent that the confidence interval associated with the eventual result will enclose the actual concentration.

It is possible to construct many different sets of intervals that satisfy this criterion. This manual uses a modification of the method of de Man (1983). De Man calculated his confidence limits iteratively from the smallest concentrations upward. Because this manual emphasizes pathogens, the intervals have been shifted slightly upward by iterating from the largest concentrations downward.

The confidence intervals of the spreadsheet and the tables associated with this appendix may be different. The MPN Excel spreadsheet uses a normal approximation to the log (MPN) to calculate its confidence intervals. This approximation is similar to a normal approximation discussed in Haldane (1939). This approximation is less computationally intense so more appropriate for a spreadsheet than de Man's confidence intervals.

Precision, Bias, and Extreme Outcomes

The MPNs and confidence limits have been expressed to 2 significant digits. For example, the entry "400" has been rounded from a number between 395 and 405.

Numerous articles have noted a bias toward over-estimation of microbial concentrations by the MPN. Garthright (1993) has shown, however, that there is no appreciable bias when the concentrations are expressed as logarithms, the customary units used for regressions and for combining outcomes. Therefore, these MPNs have not been adjusted for bias.

The outcome with all positive tubes in each dilution gives no upper bound on the concentration. The tables in this appendix list the MPN for this outcome as greater than the highest MPN for an outcome with at least one negative tube. Similarly, the outcome with all negative tubes is listed as less than the lowest MPN for an outcome with at least one positive tube.

Cautionary Notes

Improbable Outcomes

Several potential problems may cause improbable outcomes. For example, there may be interference at low dilutions or selecting too few colonies at low dilutions for a confirmation test may overlook the target microbe. If the problem is believed limited to the low dilutions, then using only the high dilutions with positive tubes might be more reliable. If the cause of the problem is unknown, then the estimate may be unreliable.

When excluding improbable outcomes, de Man's (1983) preferred degree of improbability was adopted. The outcomes included are among the 99.985 percent most likely outcomes if their own MPNs were the actual bacterial concentrations. Therefore, among 10 different outcomes, all will be found in these tables at least 99 percent of the time.

Inconclusive Tubes

In special cases where tubes cannot be judged either positive or negative (e.g., plates overgrown by competing microflora at low dilutions), these tubes should be excluded from the outcome. The resulting outcome may have different numbers of tubes than any of the tables in this appendix. Its MPN can be solved by computer algorithms or estimated by Thomas's Rule below. Haldane's method can find the confidence limits as described below Thomas's rule.

Using Tables

Selecting Three Dilutions for Table Reference

An MPN can be computed for any positive number of tubes at any positive number of dilutions, but often serial dilutions use three or more dilutions and a decimal series (Each dilution has one tenth as much of the original sample as the previous dilution.) The tables in this appendix require reducing an outcome to three of its decimal dilutions. This procedure for selecting three dilutions was developed for the designs (numbers of tubes per dilution and ratio of dilutions) in these tables. They all have decimal dilutions and a fairly small number of tubes per dilution. For other designs, other procedures may be needed. When the MPN model holds, the three decimal dilutions are chosen to give a good approximation to the MPN of the entire outcome. Otherwise, the reduction may remove interference (possible from another species of microbe or a toxic substance) that can be diluted out. The remainder of this section tells how to select the three dilutions.

First, remove the highest dilution (smallest sample volume) if it and the next lower dilution have all negative tubes. As long as this condition holds and at least four dilutions remain, continue removing these dilutions.

Next, if only three dilutions remain, use them as illustrated in example A. In each example there are five tubes in each dilution. In example A, removing the two highest dilutions (0.001 and 0.01 grams) leaves three dilutions.

If more than three dilutions remain, then find the highest dilution with all positive tubes. There are three cases. In the first case, the highest dilution with all positive tubes is within the three highest remaining dilutions. Then use the three highest remaining dilutions. In example B, the first step removes the highest dilution (0.001 grams.) Since the highest dilution with all positive tubes (1 gram) is within the three highest remaining dilutions, (1, 0.1, and 0.01 grams,) use them. In example C, the highest dilution with all positive tubes (0.01 g) is within the three highest remaining dilutions (0.1, 0.01, and 0.001.)

In the second case, the highest dilution with all positive tubes is not within the three highest remaining dilutions. Then select the next two higher dilutions than the highest dilution with all positive tubes. Assign the sum of the positive tubes of any still higher dilutions to the third higher dilution. In example D, the highest dilution with all positive tubes has 1 gram. Select the two dilutions immediately higher which have 0.1 and 0.01 grams. There is only one higher dilution whose positive tubes are assigned to

form the third dilution with 0.001 grams.

In the third case, there is no dilution with all positive tubes. Then select the two lowest dilutions. Assign the sum of the positive tubes of any higher dilutions to the third dilution. In example E no dilution has all positive tubes. The two lowest dilutions have 10 and 1 grams. The sum of the positive in the dilutions with 0.1, 0.01 and 0.001 grams is assigned to form the third dilution with 0.1 grams.

If the three dilutions selected are not in the tables, then something in the serial dilution probably was unusual. This is a warning that the outcome is sufficiently improbable that the basic assumptions of the MPN may be questionable. If possible, redoing the test may be the most reliable procedure. If an MPN value is still desired, use the three highest remaining dilutions. In example F, the three highest dilutions are used. If these dilutions are not in the tables, then use the highest dilution with any positive tubes. The section entitled 'MPN for a single dilution with any positive tubes' shows how to calculate the MPN.

Table of Examples

Examples	10 g	1 g	.1 g	.01 g	.001 g	
A	4	1	0	0	0	410xx
B	5	5	1	0	0	x510x
C	4	5	4	5	1	xx451
D	4	5	4	3	1	xx431
E	4	3	0	1	1	432xx
F	4	3	3	2	1	xx321

Conversion of Table Units

The tables below apply to inocula of 0.1, 0.01, and 0.001 g. When different inocula are selected for table reference, multiply the MPN/g and confidence limits by whatever multiplier makes the inocula match the table inocula. For example, if the inocula were 0.01, 0.001, and 0.0001 with three tubes per dilution, multiplying by 10 would make these inocula match the table inocula. If the outcome were (3, 1, 0), multiply the [Table 1](#) MPN/g estimate, 43/g, by 10 to arrive at 430/g.

Bounds and approximations for a design without a table

The MPN for a serial dilutions not addressed by any table (e.g., resulting from accidental loss of some tubes) may be computed by iteration or bounded as follows.

$$\frac{\sum_W g_j}{\sum_W \left(t_j - \frac{g_j}{2} \right) m_j + \sum_Q (t_j - g_j) m_j} \leq \lambda \leq \frac{\sum_{j=1}^K g_j}{\sum_{j=1}^K (t_j - g_j) m_j}$$

Where W and Q are two disjoint sets of dilutions that together contain all the dilutions. The lower bound allows low dilutions with all positive tubes to be deleted from the bound. Blodgett (2005b) introduces these and other bounds.

The following gives an estimate of the MPN. First, select the lowest dilution that doesn't have all positive tubes. Second, select the highest dilution with at least one positive tube. Finally, select all the dilutions between them. Use only the selected dilutions in the following formula of Thomas (1942):

$$\text{MPN/g} = (\sum g_j) / (\sum t_j m_j \sum (t_j - g_j) m_j)^{(1/2)}$$

where the summation is over the selected dilutions and

$\sum g_j$ denotes the number of positive tubes in the selected dilutions,

$\sum t_j m_j$ denotes the grams of sample in all tubes in the selected dilutions,

$\sum (t_j - g_j) m_j$ denotes the grams of sample in all negative tubes in the selected dilutions.

The following examples will illustrate the application of Thomas's formula. We assume that the dilutions are 1.0, 0.1, 0.01, 0.001, and 0.0001 g.

Example (1). For outcome (5/5, 10/10, 4/10, 2/10, 0/5) use only (–, –, 4/10, 2/10, –); so $\sum t_j m_j = 10 * 0.01 + 10 * 0.001 = 0.11$. Where * means multiplication. There are 6 negative tubes at 0.01 and 8 negative tubes at 0.001, so $\sum (t_j - g_j) m_j = 6 * 0.01 + 8 * 0.001 = 0.068$. There are 6 positive tubes, so

$$\text{MPN/g} = 6 / (0.068 * 0.11)^{(1/2)} = 6 / 0.086 = 70/\text{g}$$

Example (2). For outcome (5/5, 10/10, 10/10, 0/10, 0/5) use only (–, –, 10/10, 0/10, –), so by Thomas's formula,

$$\text{MPN/g} = 10 / (0.01 * 0.11)^{(1/2)} = 10 / 0.0332 = 300/\text{g}$$

These two approximated MPNs compare well with the MPNs for (10, 4, 2) and (10, 10, 0) (i.e., 70/g and 240/g, respectively).

Approximate confidence limits for any dilution test outcome can be calculated by first estimating the standard error of $\log_{10}(\text{MPN})$ by the method of Haldane. We describe the method for 3 dilutions, but it can be shortened to 2 or extended to any positive number.

Let m_1, m_2, m_3 denote the inoculation amounts at the largest to the smallest amounts (e.g., $m_1 = 0.1$ g, $m_2 = 0.01$ g, $m_3 = 0.001$ g in these tables).

Let g_1, g_2, g_3 denote the numbers of positive tubes at the corresponding dilutions. For legibility, we denote y^x by "y**x" and "y times x" by "y*x".

Now we compute

$$T1 = \exp(-mpn*m_1), T2 = \exp(-mpn*m_2), \text{ etc.}$$

Then we compute

$$B = [g_1*m_1*m_1*T1/((T1 - 1)**2) + \dots + [g_3*m_3*m_3*T3/((T3 - 1)**2)].$$

Finally, we compute

$$\text{Standard Error of } \log_{10}(\text{mpn}) = 1/(2.303*mpn*(B**0.5))$$

Now the 95 percent confidence intervals, for example, are found at

$$\log_{10}(\text{mpn}) \pm 1.96*(\text{Standard Error}).$$

MPN for a Single Dilution with any Positive Tubes

If just one dilution has any positive tubes, then a simpler expression gives its MPN.

$$\text{MPN}/g = (1/m)*2.303*\log_{10}((\sum t_j m_j)/(\sum (t_j - g_j) m_j))$$

Where m denotes the amount of sample in each tube in the dilution with a positive tube.

Special requirements and tables included

The attached spreadsheet should be able to handle most specialized designs. Garthright and Blodgett (2003) discusses this spreadsheet. Requests for special computations and different designs will be honored as resources permit. Designs may be requested with more or less than 3 dilutions, uneven numbers of tubes, different confidence levels, etc. (Telephone 301-436-1836 or write the Division of Mathematics, FDA/CFSAN, 5100 Paint Branch Parkway, HFS-205 Rm 2D-011, College Park, MD 20740) The most-published designs, three 10-fold dilutions with 3, 5, 8, or 10 tubes at each dilution, are presented here.

References

1. Blodgett, R. J. 2005a. Serial dilution with a confirmation step. *Food Microbiology* **22**:547-552.
2. Blodgett, R. J. 2005b. Upper and lower bounds for a serial dilution test. *Journal of the AOAC international* **88** (4):1227-1230.
3. Cochran, W. G. 1950. Estimation of bacterial densities by means of the "Most

Probable Number." *Biometrics* **6**:105-116.

4. de Man, J. C. 1983. MPN tables, corrected. *Eur. J. Appl. Biotechnol.* **17**:301-305.
5. Eisenhart, C., and P. W. Wilson. 1943. Statistical methods and control in bacteriology. *Bacteriol. Rev.* **7**:57-137.
6. Garthright, W. E. and Blodgett, R.J. 2003, FDA's preferred MPN methods for Standard, large or unusual tests, with a spreadsheet. *Food Microbiology* **20**:439-445.
7. Garthright, W. E. 1993. Bias in the logarithm of microbial density estimates from serial dilutions. *Biom. J.* **35**: 3:299-314.
8. Haldane, J.B.S. 1939. Sampling errors in the determination of bacterial or virus density by the dilution method. *J. Hygiene.* **39**:289-293.
9. Halvorson, H. O., and N. R. Ziegler. 1933. Application of statistics to problems in bacteriology. *J. Bacteriol.* **25**:101-121; 26:331-339; 26:559-567.
10. McCrady, M. H. 1915. The numerical interpretation of fermentation-tube results. *J. Infect. Dis.* **17**:183-212.
11. Peeler, J. T., G. A. Houghtby, and A. P. Rainosek. 1992. The most probable number technique, Compendium of Methods for the Microbiological Examination of Foods, 3rd Ed., 105-120.
12. Thomas, H. A. 1942. Bacterial densities from fermentation tube tests. *J. Am. Water Works Assoc.* **34**:572-576.
13. Woodward, R. L. 1957. How probable is the most probable number? *J. Am. Water Works Assoc.* **49**:1060-1068.

Tables

- [Table 1. 3 tubes each](#)
- [Table 2. 5 tubes each](#)
- [Table 3. 10 tubes each](#)
- [Table 4. 8 tubes each](#)
- [Table 5. 10 tubes at 10 ml inocula](#) (Added July 2003)

Table 1

For 3 tubes each at 0.1, 0.01, and 0.001 g inocula, the MPNs per gram and 95 percent confidence intervals.

Pos. Tubes			MPN/g	Conf. lim.		Pos. tubes			MPN/g	Conf. lim.	
0.10	0.01	0.001		Low	High	0.10	0.01	0.001		Low	High
0	0	0	<3.0	–	9.5	2	2	0	21	4.5	42

Pos. Tubes				Conf. lim.		Pos. tubes				Conf. lim.	
0.10	0.01	0.001	MPN/g	Low	High	0.10	0.01	0.001	MPN/g	Low	High
0	0	1	3.0	0.15	9.6	2	2	1	28	8.7	94
0	1	0	3.0	0.15	11	2	2	2	35	8.7	94
0	1	1	6.1	1.2	18	2	3	0	29	8.7	94
0	2	0	6.2	1.2	18	2	3	1	36	8.7	94
0	3	0	9.4	3.6	38	3	0	0	23	4.6	94
1	0	0	3.6	0.17	18	3	0	1	38	8.7	110
1	0	1	7.2	1.3	18	3	0	2	64	17	180
1	0	2	11	3.6	38	3	1	0	43	9	180
1	1	0	7.4	1.3	20	3	1	1	75	17	200
1	1	1	11	3.6	38	3	1	2	120	37	420
1	2	0	11	3.6	42	3	1	3	160	40	420
1	2	1	15	4.5	42	3	2	0	93	18	420
1	3	0	16	4.5	42	3	2	1	150	37	420
2	0	0	9.2	1.4	38	3	2	2	210	40	430
2	0	1	14	3.6	42	3	2	3	290	90	1,000
2	0	2	20	4.5	42	3	3	0	240	42	1,000
2	1	0	15	3.7	42	3	3	1	460	90	2,000
2	1	1	20	4.5	42	3	3	2	1100	180	4,100
2	1	2	27	8.7	94	3	3	3	>1100	420	–

Table 2 For 5 tubes each at 0.1, 0.01, and 0.001 g inocula, the MPNs and 95 percent confidence intervals.											
Pos. Tubes			MPN/g	Conf. lim.		Pos. tubes			MPN/g	Conf. lim.	
0.1	0.01	0.001		Low	High	0.1	0.01	0.001		Low	High

Pos. Tubes			Conf. lim.			Pos. tubes			Conf. lim.		
0.1	0.01	0.001	MPN/g	Low	High	0.1	0.01	0.001	MPN/g	Low	High
0	0	0	<1.8	–	6.8	4	0	2	21	6.8	40
0	0	1	1.8	0.09	6.8	4	0	3	25	9.8	70
0	1	0	1.8	0.09	6.9	4	1	0	17	6	40
0	1	1	3.6	0.7	10	4	1	1	21	6.8	42
0	2	0	3.7	0.7	10	4	1	2	26	9.8	70
0	2	1	5.5	1.8	15	4	1	3	31	10	70
0	3	0	5.6	1.8	15	4	2	0	22	6.8	50
1	0	0	2	0.1	10	4	2	1	26	9.8	70
1	0	1	4	0.7	10	4	2	2	32	10	70
1	0	2	6	1.8	15	4	2	3	38	14	100
1	1	0	4	0.7	12	4	3	0	27	9.9	70
1	1	1	6.1	1.8	15	4	3	1	33	10	70
1	1	2	8.1	3.4	22	4	3	2	39	14	100
1	2	0	6.1	1.8	15	4	4	0	34	14	100
1	2	1	8.2	3.4	22	4	4	1	40	14	100
1	3	0	8.3	3.4	22	4	4	2	47	15	120
1	3	1	10	3.5	22	4	5	0	41	14	100
1	4	0	11	3.5	22	4	5	1	48	15	120
2	0	0	4.5	0.79	15	5	0	0	23	6.8	70
2	0	1	6.8	1.8	15	5	0	1	31	10	70
2	0	2	9.1	3.4	22	5	0	2	43	14	100
2	1	0	6.8	1.8	17	5	0	3	58	22	150
2	1	1	9.2	3.4	22	5	1	0	33	10	100

Pos. Tubes			Conf. lim.			Pos. tubes			Conf. lim.		
0.1	0.01	0.001	MPN/g	Low	High	0.1	0.01	0.001	MPN/g	Low	High
2	1	2	12	4.1	26	5	1	1	46	14	120
2	2	0	9.3	3.4	22	5	1	2	63	22	150
2	2	1	12	4.1	26	5	1	3	84	34	220
2	2	2	14	5.9	36	5	2	0	49	15	150
2	3	0	12	4.1	26	5	2	1	70	22	170
2	3	1	14	5.9	36	5	2	2	94	34	230
2	4	0	15	5.9	36	5	2	3	120	36	250
3	0	0	7.8	2.1	22	5	2	4	150	58	400
3	0	1	11	3.5	23	5	3	0	79	22	220
3	0	2	13	5.6	35	5	3	1	110	34	250
3	1	0	11	3.5	26	5	3	2	140	52	400
3	1	1	14	5.6	36	5	3	3	180	70	400
3	1	2	17	6	36	5	3	4	210	70	400
3	2	0	14	5.7	36	5	4	0	130	36	400
3	2	1	17	6.8	40	5	4	1	170	58	400
3	2	2	20	6.8	40	5	4	2	220	70	440
3	3	0	17	6.8	40	5	4	3	280	100	710
3	3	1	21	6.8	40	5	4	4	350	100	710
3	3	2	24	9.8	70	5	4	5	430	150	1,100
3	4	0	21	6.8	40	5	5	0	240	70	710
3	4	1	24	9.8	70	5	5	1	350	100	1100
3	5	0	25	9.8	70	5	5	2	540	150	1700
4	0	0	13	4.1	35	5	5	3	920	220	2600

Pos. Tubes				Conf. lim.		Pos. tubes				Conf. lim.	
0.1	0.01	0.001	MPN/g	Low	High	0.1	0.01	0.001	MPN/g	Low	High
4	0	1	17	5.9	36	5	5	4	1600	400	4600
						5	5	5	>1600	700	–

Table 3

For 10 tubes at each of 0.1, 0.01, and 0.001 g inocula, the MPNs and 95 percent confidence intervals.

Pos. tubes			MPN/g	Conf. lim.		Pos. tubes			MPN/g	Conf. lim.	
0.1	0.01	0.001		Low	High	0.1	0.01	0.001		Low	High
0	0	0	<.90	–	3.1	8	2	0	17	7.7	34
0	0	1	0.9	0.04	3.1	8	2	1	19	9	34
0	0	2	1.8	0.33	5.1	8	2	2	21	10	39
0	1	0	0.9	0.04	3.6	8	2	3	23	11	44
0	1	1	1.8	0.33	5.1	8	3	0	19	9	34
0	2	0	1.8	0.33	5.1	8	3	1	21	10	39
0	2	1	2.7	0.8	7.2	8	3	2	24	11	44
0	3	0	2.7	0.8	7.2	8	3	3	26	12	50
1	0	0	0.94	0.05	5.1	8	4	0	22	10	39
1	0	1	1.9	0.33	5.1	8	4	1	24	11	44
1	0	2	2.8	0.8	7.2	8	4	2	26	12	50
1	1	0	1.9	0.33	5.7	8	4	3	29	14	58
1	1	1	2.9	0.8	7.2	8	5	0	24	11	44
1	1	2	3.8	1.4	9	8	5	1	27	12	50
1	2	0	2.9	0.8	7.2	8	5	2	29	14	58
1	2	1	3.8	1.4	9	8	5	3	32	15	62
1	3	0	3.8	1.4	9	8	6	0	27	12	50

Pos. tubes			Conf. lim.			Pos. tubes			Conf. lim.		
0.1	0.01	0.001	MPN/g	Low	High	0.1	0.01	0.001	MPN/g	Low	High
1	3	1	4.8	2.1	11	8	6	1	30	14	58
1	4	0	4.8	2.1	11	8	6	2	33	15	62
2	0	0	2	0.37	7.2	8	7	0	30	14	58
2	0	1	3	0.81	7.3	8	7	1	33	17	73
2	0	2	4	1.4	9	8	7	2	36	17	74
2	1	0	3	0.82	7.8	8	8	0	34	17	73
2	1	1	4	1.4	9	8	8	1	37	17	74
2	1	2	5	2.1	11	9	0	0	17	7.5	31
2	2	0	4	1.4	9.1	9	0	1	19	9	34
2	2	1	5	2.1	11	9	0	2	22	10	39
2	2	2	6.1	3	14	9	0	3	24	11	44
2	3	0	5.1	2.1	11	9	1	0	19	9	39
2	3	1	6.1	3	14	9	1	1	22	10	40
2	4	0	6.1	3	14	9	1	2	25	11	44
2	4	1	7.2	3.1	15	9	1	3	28	14	58
2	5	0	7.2	3.1	15	9	1	4	31	14	58
3	0	0	3.2	0.9	9	9	2	0	22	10	44
3	0	1	4.2	1.4	9.1	9	2	1	25	11	46
3	0	2	5.3	2.1	11	9	2	2	28	14	58
3	1	0	4.2	1.4	10	9	2	3	32	14	58
3	1	1	5.3	2.1	11	9	2	4	35	17	73
3	1	2	6.4	3	14	9	3	0	25	12	50
3	2	0	5.3	2.1	12	9	3	1	29	14	58

Pos. tubes			MPN/g	Conf. lim.		Pos. tubes			MPN/g	Conf. lim.	
0.1	0.01	0.001		Low	High	0.1	0.01	0.001		Low	High
3	2	1	6.4	3	14	9	3	2	32	15	62
3	2	2	7.5	3.1	15	9	3	3	36	17	74
3	3	0	6.5	3	14	9	3	4	40	20	91
3	3	1	7.6	3.1	15	9	4	0	29	14	58
3	3	2	8.7	3.6	17	9	4	1	33	15	62
3	4	0	7.6	3.1	15	9	4	2	37	17	74
3	4	1	8.7	3.6	17	9	4	3	41	20	91
3	5	0	8.8	3.6	17	9	4	4	45	20	91
4	0	0	4.5	1.6	11	9	5	0	33	17	73
4	0	1	5.6	2.2	12	9	5	1	37	17	74
4	0	2	6.8	3	14	9	5	2	42	20	91
4	1	0	5.6	2.2	12	9	5	3	46	20	91
4	1	1	6.8	3	14	9	5	4	51	25	120
4	1	2	8	3.6	17	9	6	0	38	17	74
4	2	0	6.8	3	15	9	6	1	43	20	91
4	2	1	8	3.6	17	9	6	2	47	21	100
4	2	2	9.2	3.7	17	9	6	3	53	25	120
4	3	0	8.1	3.6	17	9	7	0	44	20	91
4	3	1	9.3	4.5	18	9	7	1	49	21	100
4	3	2	10	5	20	9	7	2	54	25	120
4	4	0	9.3	4.5	18	9	7	3	60	26	120
4	4	1	11	5	20	9	8	0	50	25	120
4	5	0	11	5	20	9	8	1	55	25	120

Pos. tubes			Conf. lim.			Pos. tubes			Conf. lim.		
0.1	0.01	0.001	MPN/g	Low	High	0.1	0.01	0.001	MPN/g	Low	High
4	5	1	12	5.6	22	9	8	2	61	26	120
4	6	0	12	5.6	22	9	8	3	68	30	140
5	0	0	6	2.5	14	9	9	0	57	25	120
5	0	1	7.2	3.1	15	9	9	1	63	30	140
5	0	2	8.5	3.6	17	9	9	2	70	30	140
5	0	3	9.8	4.5	18	10	0	0	23	11	44
5	1	0	7.3	3.1	15	10	0	1	27	12	50
5	1	1	8.5	3.6	17	10	0	2	31	14	58
5	1	2	9.8	4.5	18	10	0	3	37	17	73
5	1	3	11	5	21	10	1	0	27	12	57
5	2	0	8.6	3.6	17	10	1	1	32	14	61
5	2	1	9.9	4.5	18	10	1	2	38	17	74
5	2	2	11	5	21	10	1	3	44	20	91
5	3	0	10	4.5	18	10	1	4	52	25	120
5	3	1	11	5	21	10	2	0	33	15	73
5	3	2	13	5.6	23	10	2	1	39	17	79
5	4	0	11	5	21	10	2	2	46	20	91
5	4	1	13	5.6	23	10	2	3	54	25	120
5	4	2	14	7	26	10	2	4	63	30	140
5	5	0	13	6.3	25	10	3	0	40	17	91
5	5	1	14	7	26	10	3	1	47	20	100
5	6	0	14	7	26	10	3	2	56	25	120
6	0	0	7.8	3.1	17	10	3	3	66	30	140

Pos. tubes			Conf. lim.			Pos. tubes			Conf. lim.		
0.1	0.01	0.001	MPN/g	Low	High	0.1	0.01	0.001	MPN/g	Low	High
6	0	1	9.2	3.6	17	10	3	4	77	34	150
6	0	2	11	5	20	10	3	5	89	39	180
6	0	3	12	5.6	22	10	4	0	49	21	120
6	1	0	9.2	3.7	18	10	4	1	59	25	120
6	1	1	11	5	21	10	4	2	70	30	150
6	1	2	12	5.6	22	10	4	3	82	38	180
6	1	3	14	7	26	10	4	4	94	44	180
6	2	0	11	5	21	10	4	5	110	50	210
6	2	1	12	5.6	22	10	5	0	62	26	140
6	2	2	14	7	26	10	5	1	74	30	150
6	2	3	15	7.4	30	10	5	2	87	38	180
6	3	0	12	5.6	23	10	5	3	100	44	180
6	3	1	14	7	26	10	5	4	110	50	210
6	3	2	15	7.4	30	10	5	5	130	57	220
6	4	0	14	7	26	10	5	6	140	70	280
6	4	1	15	7.4	30	10	6	0	79	34	180
6	4	2	17	9	34	10	6	1	94	39	180
6	5	0	16	7.4	30	10	6	2	110	50	210
6	5	1	17	9	34	10	6	3	120	57	220
6	5	2	19	9	34	10	6	4	140	70	280
6	6	0	17	9	34	10	6	5	160	74	280
6	6	1	19	9	34	10	6	6	180	91	350
6	7	0	19	9	34	10	7	0	100	44	210

Pos. tubes			Conf. lim.			Pos. tubes			Conf. lim.		
0.1	0.01	0.001	MPN/g	Low	High	0.1	0.01	0.001	MPN/g	Low	High
7	0	0	10	4.5	20	10	7	1	120	50	220
7	0	1	12	5	21	10	7	2	140	61	280
7	0	2	13	6.3	25	10	7	3	150	73	280
7	0	3	15	7.2	28	10	7	4	170	91	350
7	1	0	12	5	22	10	7	5	190	91	350
7	1	1	13	6.3	25	10	7	6	220	100	380
7	1	2	15	7.2	28	10	7	7	240	110	480
7	1	3	17	7.7	31	10	8	0	130	60	250
7	2	0	13	6.4	26	10	8	1	150	70	280
7	2	1	15	7.2	28	10	8	2	170	80	350
7	2	2	17	7.7	31	10	8	3	200	90	350
7	2	3	19	9	34	10	8	4	220	100	380
7	3	0	15	7.2	30	10	8	5	250	120	480
7	3	1	17	9	34	10	8	6	280	120	480
7	3	2	19	9	34	10	8	7	310	150	620
7	3	3	21	10	39	10	8	8	350	150	620
7	4	0	17	9	34	10	9	0	170	74	310
7	4	1	19	9	34	10	9	1	200	91	380
7	4	2	21	10	39	10	9	2	230	100	480
7	4	3	23	11	44	10	9	3	260	120	480
7	5	0	19	9	34	10	9	4	300	140	620
7	5	1	21	10	39	10	9	5	350	150	630
7	5	2	23	11	44	10	9	6	400	180	820

Pos. tubes				Conf. lim.		Pos. tubes				Conf. lim.	
0.1	0.01	0.001	MPN/g	Low	High	0.1	0.01	0.001	MPN/g	Low	High
7	6	0	21	10	39	10	9	7	460	210	970
7	6	1	23	11	44	10	9	8	530	210	970
7	6	2	25	12	46	10	9	9	610	280	1300
7	7	0	23	11	44	10	10	0	240	110	480
7	7	1	26	12	50	10	10	1	290	120	620
8	0	0	13	5.6	25	10	10	2	350	150	820
8	0	1	15	7	26	10	10	3	430	180	970
8	0	2	17	7.5	30	10	10	4	540	210	1300
8	0	3	19	9	34	10	10	5	700	280	1500
8	1	0	15	7.1	28	10	10	6	920	350	1900
8	1	1	17	7.7	31	10	10	7	1200	480	2400
8	1	2	19	9	34	10	10	8	1600	620	3400
8	1	3	21	10	39	10	10	9	2300	810	5300
						10	10	10	>2300	1300	–

Table 4. For 8 tubes each at 0.1, 0.01, and 0.001 g inocula, the MPNs per gram and 95 percent confidence intervals.

Table 4 For 8 tubes at each of 0.1, 0.01, and 0.001 g inocula, the MPNs and 95 percent confidence intervals.					
Pos. tubes			MPN/g	conf. lim.	
0.10	0.01	0.001		Low	High
0	0	0	<1.2	–	4.3
0	0	1	1.1	.057	4.3
0	0	2	2.3	.42	6.7

0.10	Pos. tubes		MPN/g	conf. lim.	
	0.01	0.001		Low	High
0	1	0	1.1	.058	4.4
0	1	1	2.3	.42	6.7
0	2	0	2.3	.42	6.7
0	2	1	3.4	1.0	9.1
0	3	0	3.4	1.0	9.1
1	0	0	1.2	.064	6.7
1	0	1	2.4	.42	6.8
1	0	2	3.6	1.0	9.1
1	1	0	2.4	.42	7.3
1	1	1	3.6	1.0	9.1
1	1	2	4.8	1.8	12
1	2	0	3.6	1.0	9.1
1	2	1	4.9	1.8	12
1	3	0	4.9	1.8	12
1	3	1	6.1	2.8	15
1	4	0	6.2	2.8	15
2	0	0	2.6	.47	9.1
2	0	1	3.8	1.0	9.1
2	0	2	5.1	1.8	12
2	1	0	3.9	1.0	9.9
2	1	1	5.2	1.8	12
2	1	2	6.5	2.8	15
2	2	0	5.2	1.8	12

0.10	Pos. tubes		MPN/g	conf. lim.	
	0.01	0.001		Low	High
2	2	1	6.5	2.8	15
2	2	2	7.9	3.3	18
2	3	0	6.6	2.8	15
2	3	1	7.9	3.3	18
2	4	0	8.0	3.3	18
2	5	0	9.4	4.3	19
3	0	0	4.1	1.2	12
3	0	1	5.5	1.9	12
3	0	2	6.9	2.8	15
3	1	0	5.6	1.9	13
3	1	1	7.0	2.8	15
3	1	2	8.4	4.0	18
3	2	0	7.0	2.9	15
3	2	1	8.5	4.0	18
3	2	2	10	4.3	19
3	3	0	8.6	4.0	18
3	3	1	10	4.3	19
3	3	2	12	5.2	24
3	4	0	10	4.3	19
3	4	1	12	5.2	24
3	5	0	12	5.2	24
4	0	0	6.0	2.1	15
4	0	1	7.5	2.9	15

0.10	Pos. tubes		MPN/g	conf. lim.	
	0.01	0.001		Low	High
4	0	2	9.1	4.1	18
4	1	0	7.6	2.9	18
4	1	1	9.2	4.1	19
4	1	2	11	4.3	22
4	2	0	9.3	4.1	19
4	2	1	11	4.3	22
4	2	2	13	5.7	24
4	3	0	11	4.3	22
4	3	1	13	5.7	24
4	3	2	14	6.6	28
4	4	0	13	5.7	24
4	4	1	15	6.6	29
4	5	0	15	6.6	29
4	5	1	16	7.2	33
4	6	0	17	7.2	33
5	0	0	8.3	3.3	18
5	0	1	10	4.3	19
5	0	2	12	5.2	24
5	0	3	14	6.6	29
5	1	0	10	4.3	22
5	1	1	12	5.2	24
5	1	2	14	6.6	29
5	1	3	16	6.7	32

Pos. tubes			conf. lim.		
0.10	0.01	0.001	MPN/g	Low	High
5	2	0	12	5.3	24
5	2	1	14	6.6	29
5	2	2	16	7.2	33
5	2	3	18	7.2	33
5	3	0	14	6.6	29
5	3	1	16	7.2	33
5	3	2	18	7.2	33
5	4	0	16	7.2	33
5	4	1	18	7.6	33
5	4	2	21	9.0	39
5	5	0	19	7.6	33
5	5	1	21	9.0	39
5	6	0	21	9.0	39
6	0	0	11	4.3	24
6	0	1	13	5.7	25
6	0	2	16	6.6	32
6	0	3	18	7.2	33
6	1	0	14	5.8	29
6	1	1	16	6.6	32
6	1	2	18	7.2	33
6	1	3	21	9.0	39
6	2	0	16	6.7	33
6	2	1	18	7.4	33

0.10	Pos. tubes		MPN/g	conf. lim.	
	0.01	0.001		Low	High
6	2	2	21	9.0	39
6	2	3	23	11	50
6	3	0	19	7.6	35
6	3	1	21	9.0	39
6	3	2	24	11	50
6	3	3	27	12	53
6	4	0	21	9.0	40
6	4	1	24	11	50
6	4	2	27	12	53
6	6	1	31	13	69
6	7	0	31	13	69
7	0	0	16	6.6	32
7	0	1	18	7.2	33
7	0	2	21	9.0	40
7	0	3	25	11	50
7	1	0	19	7.9	39
7	1	1	22	9.0	40
7	1	2	25	11	50
7	1	3	29	12	54
7	2	0	22	9.0	45
7	2	1	25	11	51
7	2	2	29	13	68
7	2	3	33	13	69

0.10	Pos. tubes		MPN/g	conf. lim.	
	0.01	0.001		Low	High
7	3	0	26	11	53
7	3	1	30	13	68
7	3	2	34	13	69
7	3	3	39	17	91
7	4	0	30	13	69
7	4	1	35	13	69
7	4	2	39	17	91
7	4	3	45	18	101
7	5	0	36	14	75
7	5	1	40	17	91
7	5	2	46	18	101
7	5	3	52	21	117
7	6	0	41	17	91
7	6	1	47	21	117
7	6	2	53	21	117
7	6	3	59	24	146
7	7	0	48	21	117
7	7	1	55	21	117
7	7	2	61	24	146
7	8	0	56	21	119
8	0	0	23	9.7	50
8	0	1	28	12	54
8	0	2	34	13	69

0.10	Pos. tubes		MPN/g	conf. lim.	
	0.01	0.001		Low	High
8	0	3	41	17	91
8	1	0	29	12	68
8	1	1	35	13	75
8	1	2	43	17	91
8	1	3	52	21	120
8	1	4	63	28	150
8	2	0	36	14	91
8	2	1	45	17	100
8	2	2	55	21	120
8	2	3	67	28	150
8	2	4	81	32	190
8	3	0	47	18	120
8	3	1	58	21	150
8	3	2	72	28	150
8	3	3	87	39	190
8	3	4	102	39	190
8	3	5	118	50	240
8	4	0	62	24	150
8	4	1	77	28	190
8	4	2	94	39	190
8	4	3	110	44	220
8	4	4	130	53	250
8	4	5	150	68	280

0.10	Pos. tubes		MPN/g	conf. lim.	
	0.01	0.001		Low	High
8	5	0	84	32	190
8	5	1	100	39	220
8	5	2	120	50	250
8	5	3	140	67	280
8	5	4	170	74	340
8	5	5	190	74	340
8	5	6	210	90	400
8	6	0	110	45	240
8	6	1	140	53	280
8	6	2	160	68	340
8	6	3	190	74	340
8	6	4	220	90	400
8	6	5	250	120	490
8	6	6	290	120	520
8	7	0	160	68	340
8	7	1	190	74	400
8	7	2	230	90	490
8	7	3	270	116	520
8	7	4	310	150	710
8	7	5	370	150	720
8	7	6	430	190	1000
8	7	7	510	190	1000
8	8	0	240	99	490

Pos. tubes			MPN/g	conf. lim.	
0.10	0.01	0.001		Low	High
8	8	1	300	120	710
8	8	2	380	150	1000
8	8	3	510.0	190	1200
8	8	4	700	240	1700
8	8	5	980	340	2200
8	8	6	1400	490	3100
8	8	7	2100	710	5100
8	8	8	>2100	1000	–

(Added July 2003)

Pos. tubes	MPN/100ml	Conf. lim.	
		Low	High
0	<1.1	–	3.3
1	1.1	.05	5.9
2	2.2	.37	8.1
3	3.6	.91	9.7
4	5.1	1.6	13
5	6.9	2.5	15
6	9.2	3.3	19
7	12	4.8	24
8	16	5.9	33
9	23	8.1	53

Conf. lim.

Pos. tubes	MPN/100ml	Low	High
10	>23	12	–

(Added July 2003)

Download (</downloads/Food/FoodScienceResearch/UCM164329.zip>) an Excel spreadsheet to calculate values (zip file).

Note: The confidence intervals of the spreadsheet and the tables associated with this appendix may be different. The MPN Excel spreadsheet uses a normal approximation to the log (MPN) to calculate its confidence intervals. This approximation is similar to a normal approximation discussed in Haldane (1939). This approximation is less computationally intense so more appropriate for a spreadsheet than de Man's confidence intervals.

Hypertext Source: BAM 8th Edition, Modified from Revision A CD ROM version 1998 on 6/21/2000.

More in Laboratory Methods
[\(/Food/FoodScienceResearch/LaboratoryMethods/default.htm\)](/Food/FoodScienceResearch/LaboratoryMethods/default.htm)

[CFSAN Laboratory Quality Assurance Manual \(/Food/FoodScienceResearch/LaboratoryMethods/ucm2006948.htm\)](/Food/FoodScienceResearch/LaboratoryMethods/ucm2006948.htm)

[Microbiological Methods & Bacteriological Analytical Manual \(BAM\) \(/Food/FoodScienceResearch/LaboratoryMethods/ucm114664.htm\)](/Food/FoodScienceResearch/LaboratoryMethods/ucm114664.htm)

[Drug & Chemical Residues Methods \(/Food/FoodScienceResearch/LaboratoryMethods/ucm2006950.htm\)](/Food/FoodScienceResearch/LaboratoryMethods/ucm2006950.htm)

[Elemental Analysis Manual \(EAM\) for Food and Related Products \(/Food/FoodScienceResearch/LaboratoryMethods/ucm2006954.htm\)](/Food/FoodScienceResearch/LaboratoryMethods/ucm2006954.htm)

[Macroanalytical Procedures Manual \(MPM\) \(/Food/FoodScienceResearch/LaboratoryMethods/ucm2006953.htm\)](/Food/FoodScienceResearch/LaboratoryMethods/ucm2006953.htm)

[Pesticide Analytical Manual \(PAM\) \(/Food/FoodScienceResearch/LaboratoryMethods/ucm2006955.htm\)](/Food/FoodScienceResearch/LaboratoryMethods/ucm2006955.htm)